

FVS Kit Supplies

- Acrylic tube
 - Clear is best for teaching but white PVC could also work
 - 6" diameter, approx. 6-8" tall, can purchase online or at hardware store
 - Drill 4 holes at 12, 3, 6 and 9 o'clock.
 - Can do several sets of four holes at different depths along the tube to practice superficial vs deep suturing
 - Optional: non-stick surface of some type so the tube doesn't slide. Rubber hose along the base of the tube, drawer liner, etc.
- Alligator clips
 - Must fit through the holes drilled in the acrylic tube, claw facing inside
 - Attach to the outer side of the acrylic tube to maintain tension on the foam sheet. Elastic string, rubber bands, etc.
- Foam sheets
 - ¼" thickness, multiple sizes sold at any craft store, online or target, etc.
 - Cut into 4" x 4" squares
 - Stamp or mark with markers per the clockface pattern
- Stamps
 - Two stamps, can be custom ordered from Etsy for \$20-30
 - Use dimensions from clockface pattern and patch pattern (see below)
 - Can also trace onto foam using a marker and a pre-measured pattern
- Graft-type material
 - Long balloons (the type used to make balloon animals)
 - Gloves (ultra-thin can present an extra challenge!)
 - Odds and ends from vascular cases (pass off extra graft before it touches the patient!)
 - Ask local Gore rep for a donation
 - 8mm PTFE for end to side, 10mm PTFE for patch
- Suture
 - 3-0 Prolene with SH needle for clockface
 - 5-0 C1 prolene for end to side and patch
 - Ask local Ethicon rep for educational donation, unused sutures from OR

Clock Face Stamp Dimensions

